

Reverse Transcription protocol

Eric Greer 9/2/09

Take 2 ug of total RNA sample in 9 λ total volume

Add 0.5 λ Dnase

Add 0.5 λ RNase OUT (I make a master mix of the two and add 1 λ)

Put @ 37° for 15 minutes

Move to 75° for 15 minutes

Quick spin

Keep 9 λ for cDNA and remove 1 λ as control for RNA (will add 9 λ dH₂O to equalize volume)

To 9 λ add 1 λ dNTP and 1 λ OligodT primers (again I add 2 λ of a master mix)

Put at 65° for 5 minutes

Put on ice

Quick spin

Add 9 λ of master mix

4 λ 25mM MgCl₂

2 λ 10X RT buffer

2 λ 0.1M DTT

1 λ RNase OUT

9 λ

incubate at 42° for 2 minutes

Add 1 λ SSII

Pipette gently up and down to mix

Put at 42° for 50 minutes

Move to 70° for 15 minutes

Quick spin

Add 80 λ dH₂O Keep @ -20 use 3.3 λ per Real Time reaction